miniland

89718 easy flow duo

Español 6
English 20
Português 33

Français 47
Deutsch 61
Italiano 76

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1. INTRODUCTION

We congratulate you on purchasing this lightweight, portable, and hands-free dual breast pump from Miniland, designed to provide a more efficient expression of breast milk. Before using this device, please carefully read all instructions in this manual. Keep the manual for future reference as it contains important information.

The features described in this user guide are subject to change without notice.

2. CONTENT

- 2 Main unit
- 2 24mm breast shield
- 2 15mm, 17mm, 19mm and 21mm adaptors
- 2 USB cable
- 2 Bra adjustment buckle
- 2 Antidust cover
- 1 Nipple measure card
- 1 Carrying bag
- 1 Instruction manual
- 1 Guarantee card

If any of the above items are missing, please, contact your distributor.

3. IMPORTANT SAFETY INFORMATION

Please read all instructions before using this product.
Keep these instructions for future reference.

Warnings identify all instructions that are important for safety. Failure to observe these instructions can lead to personal injury or damage to the product.

When using electrical products, especially when children are present, basic safety precautions must always be followed.

WARNING

To avoid fire, electric shock, or serious burns:

- The breast pump should never be left unattended when plugged into a power source.

ENGLISH

- Do not use an electrical product that has been exposed to water or other liquids including:
 - Do not use while bathing or showering.
 - Never place or drop into water or other liquids.
 - Do not run water over the breast pump.
 - If an electrical product has been exposed to water or other liquids, do not touch and if possible unplug the device from electrical outlet, turn off and contact the manufacturer.
- Be aware of the risk of USB-C-socket of the pump or battery being shortcircuited by metal objects or fluids.

CAUTION

To avoid health risks and reduce the risk of severe injury:

- This is a single-user product. Use by more than one person may present a health risk and voids the warranty.
- Do not use the breast pump while operating a moving vehicle.
- Pumping may induce labour. Do not pump until after giving birth. If you become pregnant while breastfeeding or breast pumping, consult with a licensed healthcare professional before continuing.
- If infected with Hepatitis B, Hepatitis C, or Human Immunodeficiency Virus (HIV), pumping breast milk will not reduce or remove the risk of transmitting the virus to your baby through your breast milk.
- Wash all parts that come into contact with your breast and breast milk after every use.
- Be aware that small parts may be swallowed by toddlers.
- Close supervision is necessary when the breast pump or accessories are used in the vicinity of children.
- This appliance is not intended for use by persons (including children up to 14 years) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children must not play with the appliance.
- Cleaning and maintenance must not be performed by children without supervision.

 CAUTION**To avoid fire or burns:**

- Always unplug electrical products immediately after charging.
- Do not place or store the product where it can fall into water.
- The breast pump and accessories are not heat-resistant. Keep away from heated surfaces or open flames.
- Do not use near flammable materials.

 CAUTION**To avoid health risks and reduce the risk of injury:**

- Before each use visually inspect the individual components for cracks, chips, tears, discoloration or deterioration. In the event that damage to the device is observed, please discontinue use until the parts have been replaced.
- Do not repair yourself. Do not modify the device.
- Never use a damaged device. Replace defective or worn parts.
- Use the breast pump only for its intended purpose as described in these instructions for use.
- Do not use the breast pump while sleeping or overly drowsy.
- Do not microwave or boil breast milk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving (microwaving or boiling can also change the composition of breast milk).
- Clean and sanitise all parts that come into contact with your breast and breast milk prior to first use.
- If pumping is uncomfortable or causing pain, turn the unit off, break the seal between the breast and the breast shield with your finger and remove the breast shield from your breast.

4. PRODUCT FEATURES**4.1. COMPONENTS (FIGURE 1)**

1. Breast shield 24mm
2. Nipple Tunnel
3. Silicone diaphragm
4. Linker
5. Valve

■ ENGLISH

6. Milk container
7. Pouring opening
8. Motor unit
9. USB cable
10. Bra adjustment buckle
11. USB charging port
12. Antidust cover
13. Nipple measure card

4.2. BUTTONS (FIGURE 2)

- I. Function button
- II. Start/Stop button
- III. Increase vacuum button
- IV. Decrease vacuum button

4.3. DISPLAY (FIGURE 3)

- a. Vacuum level display
- b. Session time in minutes
- c. Battery display
- d. Expression phase indicator
- e. Stimulation phase indicator
- f. Unit of running time

5. INSTRUCTIONS FOR USE

5.1. INTENDED PURPOSE

The electric breast pump is a powered breast pump to be used by lactating women to express and collect milk from their breasts. The electric breast pump is intended for a single user. This breast pump is intended to be used in household, offices, hotels and other residential type environments.

5.2. PRODUCT DESCRIPTION

The electric breast pump is a personal-use electric breast pump that includes two operational modes (Stimulation and Expression) and is capable of single pumping.

5.3. BEFORE USE

5.3.1. POWERING YOUR BREAST PUMP

Before first time use, help to charge the battery of the device. Plug the power cable into the supply socket on top of the motor unit (11) and charge the battery for two hours. Note that the battery icon on the display (c) flashes when charging.

ATTENTION: Use only an external adapter with the following specifications: <input power AC 100~240V 50/60 Hz and output power is DC 5V 1.0A>

WARNING:

- Do not expose the breast pump or battery to excessive temperatures.
- The permissible ambient temperature for charging the battery is 5°C - 35°C.
- This appliance contains batteries that are non-replaceable.
- When the battery is at end of life, the appliance shall be properly disposed of.
- The integrated battery must not be removed for disposal.
- If the battery should leak, the breast pump must be disposed of.

5.3.2. DISSASSEMBLING THE PRODUCT

Before the first use of the device, all parts in contact with the milk must be washed, disinfected and dried (see chapter 5.3.3). To do it, the device should be disassembled following the next steps (**Figure 4**):

1. Open the collection cup using the tabs at the bottom.
2. Remove the linker from the breast shield.
3. Remove the silicone diaphragm and valve.

5.3.3. CLEANING

It is important you clean the complete product before using your breast pump for the first time. To do it, please follow the next instructions.

5.3.3.1. CLEANING SILICONE AND PLASTIC PARTS

For the breast diaphragm, milk container, linker and valve, please clean them with the indicated frequency:

- Before first use, wash and sanitise the components.
- After each use, wash the components.
- Once per day after use, wash and sanitise the components.

■ ENGLISH

5.3.3.1.1. REGULAR CLEANING

- You can wash these parts using the dishwasher and soft dishwashing liquid. Do not use disinfectant, as it could damage the product.
- You can wash these parts by hand using a clean wash basin used for infant items. Do not place these parts in the sink to rinse or wash.
- It is recommended to use drinking-quality tap water or bottled water.
- Take care when cleaning the valve and silicone diaphragm, if it gets damaged, your breast pump will not function properly.

5.3.3.1.2. STERILIZATION

You can use different methods to sanitise the silicone and plastic parts of your product.

Sterilization in boiling water:

Cover the disassembled, washed parts, with plenty of water and boil at least for 10 minutes. Ensure that the parts are not in direct contact to the bottom of the pot.

Sterilization in microwave:

Use bags in the microwave in accordance with the instructions on the bags.

Notes:

- If you notice a white residue on your collection cups parts after boiling, your water may have a high mineral content. Remove residue by wiping parts with a clean towel and allow to air dry.
- Distilled water is recommended when boiling parts to prevent substantial mineral build-up over time, which may compromise your parts.

5.3.3.2. MOTOR UNIT CLEANING

Wipe with clean, damp towel. Other cleaners or abrasive cloths may degrade the display performance. Do not immerse the motor unit in water; do not run water over the motor unit.

5.3.4. DRYING

1. Allow to air dry on a clean, unused dish towel or paper towel.
2. Put the clean and dry parts in a clean storage bag or a clean environment. Do not store parts in an airtight container/bag if moist.

3. Store the dry collection cups in a clean bag/container until next use.

5.3.5. ASSEMBLY OF THE BREAST PUMP

To prevent damage to the breastpump all components must be completely dry before use.

If after assembling, you don't feel suction, please check that you have followed the next steps correctly (**FIGURE 5**):

1. Put the silicone diaphragm into the linker and press it with your hands to make the edges fasten tightly to ensure a perfect seal.
2. Install the valve into the bottom of the linker, make sure it is installed in the correct direction.

Note 1: As an improperly seated valve will prevent the breast pump from generating suction.

3. Put the assembled linker into the milk container, with the "ear" part of the linker aligned with the internal slot of the milk container and press the linker disc part with both hands to fasten the edges firmly to ensure complete sealing.
4. Install the silicone shield on the milk container, align the deflector on the silicone shield with the deflector on the milk container, and press the edges with both hands to ensure that it fits into place.
5. Turn the motor unit off, connect the motor unit into the milk collector ensure that it fits perfectly.
6. Press the silicone shield firmly against the breast. Ensure that the nipple is centered in the tunnel of the breast shield. Be careful not to distort the position of the product and the chest, otherwise the breast pump will not have suction. Place the collection cup in your bra.

Note 2: If you experience pain or discomfort, the device may not be properly aligned. Please adjust the alignment for comfort.

Note 3: Make sure the breast pump and the milk container always stay at a vertical angle. Otherwise, the milk will spill out.

7. Secure the collection cups using a nursing bra and press the pump firmly against the breast. Tighten the bra or use the adjustment buckle to ensure a secure fit.

ENGLISH

CAUTION

To prevent contamination of your milk:

- Wash hands thoroughly with soap and water before touching the motor unit, parts of the collection cups and breasts.
- Avoid touching the inside of the collection cups components.

To make sure that your breast pump works properly and safely:

- Check collection cups components for wear or damage before use. Replace if necessary.
- Always inspect all parts prior to use for cleanliness. If dirty, clean according to chapter 5.3.3.

5.4. OPERATING YOUR BREAST PUMP

5.4.1. TERMS AND DESIGNATIONS

Stimulation function	Fast sucking/pumping mode to stimulate milk flow.
Expression function	Slower sucking/pumping mode for gentle and efficient milk removal after milk has started flowing.
Maximum comfort vacuum	Highest vacuum setting where pumping still feels comfortable. Different for every mother.
Let-down	The reflex that makes the milk begin to flow from the breast.

5.4.2. MILK EXTRACTION

CAUTION

To extract the milk using the breastpump, please, follow the next steps:

1. Ensure the collection section are correctly connected and positioned.
2. Press and hold the Start/Stop button (II) to switch on the device.
3. Select the operating mode (stimulation or expression) by pressing the mode selection button (I) and then press the Start/Stop button (II) again to begin the session.

4. During the early postpartum period, the milk ducts may not yet be fully open; in that case, you can use the stimulation function (e) to stimulate the breast and promote the flow of breast milk. As soon as the milk starts to flow, you can press the mode selection button (I) to select the expression mode (d).
5. Adjust the vacuum with the Increase vacuum (III) and Decrease vacuum (IV) buttons to find a comfortable level. Once you are pumping in the Expression phase, increase the vacuum level until pumping feels slightly uncomfortable (not painful), then decrease the vacuum by 1 level with the Decrease vacuum (4) button.
6. If you need to temporarily stop pumping, briefly press the Start/Stop button (II). This allows you to readjust the position of the breast pump or attend to your surroundings. If you do not resume pumping within 3 minutes, the device will switch off automatically.
7. When you finish the pumping session completely, stop the breast pump by pressing and holding the Start/Stop button (II). The breast pump will automatically switch off after 30 minutes if there is no user interaction.

Information

- Maximum comfort vacuum is the highest vacuum setting where pumping still feels comfortable. This is different for every mother.
- Stimulation should be at a comfortable vacuum level, striving for a maximum level is not necessary.
- Reassess your maximum comfort vacuum throughout your pumping experience. It may change throughout each stage of lactation.

CAUTION

- Avoid using a vacuum setting that is too high and causes discomfort or pain. The pain, along with potential breast and nipple trauma, may decrease milk output.
- Contact your healthcare professional or breastfeeding specialist if you can express only minimal or no milk or if expression is painful.

5.5. HANDLING OF BREAST MILK

5.5.1. POURING MILK (FIGURE 6)

1. While the pump is still running, slightly lean forward and ensure all milk from the nipple tunnel enters the collection container. Turn off pump.

ENGLISH

2. Carefully remove the collection cups from your bra.
3. You can use the volume scale markings on the shield to read the collected milk volume. The scale on the milk container gives you an approximate indication of the volume during pumping.
4. Hold the milk collection container facing downwards and pour the breast milk through the pouring opening (7) at the top of the container directly into a baby bottle or a clean storage container. Do not use the collection container to store milk or keep it assembled in the refrigerator.

CAUTION

Separate and wash all parts that are exposed to breast milk immediately after use. This will help to remove breast milk residue and prevent growth of bacteria.

5.5.2. GENERAL STORAGE GUIDELINES

Storage guidelines for freshly expressed breast milk (for healthy term babies)			
Room temperature 16 to 25°C (60 to 77°F)	Refrigerator 4°C (39°F) or colder	Freezer -18°C (0°F) or colder	Breast milk thawed in the refrigerator
Up to 4 hours (under very hygienic conditions and at low room temperature, up to 6 h, although not ideal)	Up to 3–4 days (maximum 5 days if strict hygiene measures have been followed)	Up to 6 months is ideal *Up to 12 months if milk was expressed under strict hygienic conditions	At room temperature: up to 2 hours In the refrigerator: up to 24 hours Do not refreeze!

NOTICE

- Repeated freeze/thaw cycles destroy the structure of the milk. Therefore, never refreeze breast milk.

- Freeze expressed breast milk in plastic food-grade bottles or milk storage bags. Do not fill the bottles or bags more than 3/4 full to allow space for possible expansion.
- Label the bottles or bags with the date and volume of expression.
- Plastic bottles and parts become brittle when frozen and may break when dropped.
- Do not use the breast milk from bottles or components that show signs of damage.

 CAUTION

- Always check the milk temperature before feeding.
- Do not thaw or warm breast milk in a microwave or a pan of boiling water. This helps preserve important components and prevents burns.
- To preserve breast milk components, thaw the milk overnight in the refrigerator.
- Alternatively, hold the bottle or bag under warm water (max. 37°C) and use immediately, as soon as it has been thawed.

5.6. TROUBLESHOOTING

Problem	Possible cause	Solution
The timer does not count up after you pressed the Start/Pause button	The product is defective	Contact Miniland customer service
The breast pump generates no vacuum (motor not working) after you pressed the Start/Pause button	The breastpump is not turned on	Make sure that the breast pump is switched on.
	The battery is discharged	Charge the device as indicated on paragraph 5.3.1.
	The product is defective	Contact Miniland customer service

ENGLISH

Problem	Possible cause	Solution
There is low or no suction	The components, specially the valve is not clean	Make sure that all collection cups components are clean and dry and in the correct direction, especially valves.
	The breast shield is not correctly positioned	Make sure the breast shields form a complete seal around the breast.
	The assembly of the product is not correct	Follow carefully the assembly instructions indicated on paragraph 5.3.5
	Misalignment of the main silicone insert.	Make sure the funnel size is well chosen and the nipple is in the tunnel.
	The product is defective	Contact Miniland customer service
Milk leakage	The breast shield is damaged or has been incorrectly fitted.	Stop pumping immediately, switch off and wash the device, and fit the membrane correctly or replace it.

6. MAINTENANCE AND CARE

6.1. BATTERY LIFE

To preserve the battery life cycle, please follow the next recommendations:

- Keep your breast pump in a cool place.
- Recharge the battery before it completely drains. This is better for the life cycle of the battery.
- Fully charge the battery before storing your breast pump.

6.2. CARE OF THE PRODUCT

To preserve all the components and parts on good conditions for as long as possible, please:

- Clean it according to the instructions on paragraph 5.3.3.
- Place the dust cover on the product once it has been cleaned and it's completely dry
- Place it inside the dust bag when it's not in use to prevent the product from getting dirty
- Keep the product on a dry and cool environment while it's not in use.

7. TECHNICAL SPECIFICATIONS

Rated Power	1.0W
Core Functions	Massage/Lactation
Levels	9 Levels
Display type (Size)	LED (25×15MM/1.2inch)
Milk collection cup capacity	180ml/6oz
Negative Pressure Range (Massage Mode)	level 1> -30mmHg, level 9> -130mmHg
Negative Pressure Range (Lactation/ Expression Mode)	level 1≥ -45mmHg, level 9 ≥ -210mmHg
Maximum Negative Pressure Range	max. -245mmHg
Negative Pressure Pump Motor Speed	6100r/min
Noise Level	<55dBa
Battery Capacity	DC 3.7V/1000mAh lithium battery
Charging Time	150 mins
Working Time (Lactation Level 9)	130 mins
Automately turn off	Standby: 5 mins Operating 30 mins

ENGLISH

Waterproof Rating	IP22
Dimensions (L×W×H)	74×111×141mm
Weight	215g

7.1. MEANING OF SYMBOLS

Direct current

Caution

Manufacturer

Date of manufacture

Batch code

Comply with local regulations about disposal

Fragile handle with care

Keep away from rain

Keep away from sunlight.

IP22 2 Protected against solid foreign objects of 12.5 mm Ø and greater
2 Protected against vertically falling water drops when enclosure tilted up to 15°

8. PRODUCT DISPOSAL INFORMATION

8.1. DISPOSAL OF YOUR BREAST PUMP

At the end of its operating life, separate the parts of your breast pump and dispose of according to the following instructions.

8.1.1. MILK CONTAINER AND LINKER

The parts are made of plastics that are not harmful to the environment when disposed of as household waste. Recycle or dispose of according to local regulations.

8.1.2. MOTOR UNIT

- It is necessary to dispose of the product responsibly.
- Do not dispose of products showing the crossed-out wheelie bin with other household waste.
- Once the product's service life has ended, it must be disposed of by taking it to a designated collection and recycling center. You may also contact the store where the product was purchased.

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We declare under our own responsibility that the product:

easy flow duo · Miniland · 89718 · Miniland, S. A.

**to which this declaration refers conforms with the
relevant standards or other standardising documents:**

- **Electrical safety (Standard 2014/35/EU):**

- EN 60335-1:2012 + A11:2014 + A13:2017 + A1:2019 + A2:2019 + A14:2019
Household and similar electrical appliances – Safety – Part 1: General requirements
- EN 62233:2008 Measurement methods for electromagnetic fields of household appliances and similar apparatus with regard to human exposure

- **Electromagnetic Compatibility (Directive 2014/30/EU):**

- EN IEC 55014-1:2021 Electromagnetic compatibility – Emission requirements for household appliances
- EN IEC 55014-2:2021 Electromagnetic compatibility – Immunity requirements
- EN IEC-3-2:2019/A1:2021 Limits for harmonic current emissions (equipment input current ≤ 16 A per phase)
- EN 61000-3-3:2013/A2:2021 Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems

- **Regulation 1935/2004/EC on materials and articles intended to come into contact with food**

- **Regulation (EU) No.10/2011 on plastic materials and articles intended to come into contact with food**

DECLARATION OF CONFORMITY

Miniland S. A. P.Ind. La Marjal I C/ La Patronal 10. 03430 ONIL (Alicante) SPAIN
Technical support: +34 966 557 775 · www.minilandgroup.com · miniland@miniland.es

We declare under our own responsibility that the product:

easy flow duo · Miniland · 89718 · Miniland, S. A.

**to which this declaration refers conforms with the
relevant standards or other standardising documents:**

- Resolution AP(2004)5 on silicones used for food contact applications
- Restriction of Hazardous Substances (RoHS) 2.0 Directive 2011/65/EU Annex II (EU) 2015/863 as last amended by Directive (EU) 2017/2102
- Regulation (EU) 2023/1542 of the European Parliament and of the Council of 12 July 2023 concerning batteries and waste batteries, amending Directive 2008/98/EC and Regulation (EU) 2019/1020 and repealing Directive 2006/66/EC
- REACH (Regulation concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals) Regulation (EC) No 1907/2006
- Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE)

Miniland, S.A.
Chief Executive Officer
Julien Fournier

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Made in China. Keep this information.

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© Miniland, S. A. 2025
Miniland S.A. · P.Ind. La Marjal · C/ La Patronal 10