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
BATHROOM SCALE

MOD. BARANOA 309-10886

Please read the instructions carefully before use to ensure safe and satisfactory operation of this product.

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.

Model: iF2832WA

Capacity: Max 180kg/400lb

Division: $100\text{g} \leq \text{CAPACITY} < 10\text{kg}$ $d=10\text{g}/0.02\text{lb}$

$10\text{kg} \leq \text{CAPACITY} \leq 180\text{kg}$ $d=50\text{g}/0.1\text{lb}$

BRIEF INTRODUCTION

This scale adopts 8 electrodes, multi-frequency BIA and other technologies, which can accurately measure and analyse a full range of body data. Used in conjunction with the APP, it can display and record comprehensive body data, as well as analyse and draw the change trajectory, and check the body changes at any time. It's the first step to control your weight and keep your body shape.

STATEMENT

The pictures, colours, accessories, user interface and other illustrations in the manual are for illustration only. Due to product upgrades, there may be slight differences between the actual product

IMPORTANCE INSTRUCTIONS

Accurate body data may not be available for the following groups of people:

- Minors, seniors over 80 years old.
- Professional athletes.
- Pregnant women, artificial dialysis patients, people in a swollen state.

Individuals with pacemakers or other implanted medical devices should only use this product in Zero-Current Mode.

This product is not a medical device. Any measurement data provided is for reference only and should not be used as medical advice.

MODE INSTRUCTIONS

- Weight measurement mode: No need to pull up the handle. Step on the scale to complete the measurement.
- Zero-Current Mode: When enabled, the bioimpedance current is turned off. The device will only measure weight and BMI, and the display will show the “” icon. If you are pregnant or have implanted medical devices, please enable this mode.

1




Zero-Current Mode is disabled by default and can be enabled via the App.
App Home > Device Page > Zero-Current Mode.

- Offline Mode: Once you have completed the initial pairing process and taken a first measurement, you will not need to connect your phone again to save each subsequent measurement. The data will automatically upload to the App the next time you connect your phone.
- Small object weighing mode: The scale can directly weigh objects over 0.1kg. For small objects under 10kg, you need to step on the scale surface first to wake up the scale (apply a force greater than 5kg and press the scale surface) and wait until the scale has returned to zero before weighing.


2

The Small Item Weighing Mode is disabled by default and can be enabled via the App.
App Home > Device page > Small Item Weighing Mode.


The usage instructions after enabling Small Item Weighing Mode are as follows:
After installing batteries or moving the weighing scale, the display may show the scale's own weight (Nothing has been weighed, but the display shows a weight). Wait six seconds for the display to return to zero before use.

GENERAL VIEW OF THE PRODUCT - FRONT SIDE


GENERAL VIEW OF THE PRODUCT - BACK SIDE

BATTERY INSTALLATION

4*1.5V AA batteries are used for this scale.
Please replace them as per the instructions below:
Change battery: Open the battery cover, take out the used batteries, then install the new ones properly as shown below.

Note:
Please use the right battery. Remove battery if the scale leaves unused for more than 3 months.
The battery and related pictures are for reference only, please refer to the actual product.

HOW-TO-USE


Step 1: Place the scale properly and download the APP.

Search “iMoni Health” in “App Store” or “Google Play” to download the Application.
Or in another most convenient way by scanning the QR code above.

Step 2: Register and bind the device

Turn on the Bluetooth on your phone, open the APP, and follow the guidelines to register and bind the device:

- If you have purchased a WiFi version, follow the App's instructions during the initial setup process to connect to Wi-Fi.
If you missed the initial connection or need to reconfigure your Wi-Fi, please follow the steps below:

App Home > Devices page > Wi-Fi Configuration.

Warm reminder:

Personal body information is only used for body composition analysis, not for other purposes. To ensure the accuracy of the results, please enter accurate personal information.


Step 3: Measure

• Body Composition Analysis Measurements


1. Place the scale on a hard, flat surface and open the App.
2. Gently step on the scale with bare feet and place your feet on the 4 electrodes on the left and right side of the scale, then lift up the handle. Please hold the handle with both hands as shown in the picture, adjust the height of both hands and straighten your arms (it is best if your arms don't touch your waist).


Step on the scale with bare feet and place your feet on the 4 electrodes


Hold the two front electrodes with thumbs


Hold the two back electrodes with the rest of the fingers

• Wrong Measuring Postures

The ten fingers are away from the electrodes


The two arms are not straight or stretched

Squatting down to measure

The knees are bent or not straightened

Measuring with socks

The two feet fail to step on the 4 electrodes evenly

3. After the weight value displayed in the window has stabilised and is blinking, the body composition analysis measurement begins.

Body weight measurement

Body fat, body muscle, etc. measuring

Measurement complete

4. Check the measurement reports with the APP

View the measurement reports in the APP and save the historical data.

Tips:

- The whole process will take about 20 seconds.
- If the APP is not connected for the first time of measuring, only weight can be measured, but no body fat. Please connect the APP before using.
- App is kept updating, so it may differ from showed on manual.
For the final and accurate information, please check on App.

MULTIPLE USERS SUPPORTED

Multiple users binding in one scale is supportable. A maximum of 15 groups of historical data for each user can be stored. Data will be automatically synchronized to the smart device when the APP is running.

Note: Final user amount depends on the product.

DISPLAY ICONS


(Note: Display differs from one to another, so below picture is for indicative only.)

	Fat		Muscle
	Wifi		Bluetooth
	Zero-Current Mode		%
	Weight Unit		Weight Unit
	Weight Unit		

	Measuring
	Wrong Measurement Indication
	Low Battery Indication
	Over-load Indication
	Negative Weight Alert
	Re-Test Indication

ABOUT THE UNIT

- It is suggested to set the unit on App.
- The scale has a unit conversion button. Open the battery cover at the bottom of the scale, in the position as shown below, and press this button to select the unit you need.
- The factory unit setting differs from different sales destination. Please refer to the actual product.

ADVICE FOR USE & CARE

1. To ensure the accuracy of the data, it is recommended that each measurement be taken at the same time of day, with the recommended time of measurement being 2 hours after a meal.
 - * Body weight changes with eating and drinking, and body fat percentage = fat weight/weight, so it is normal for fat percentage to fluctuate.
2. Please do not make measurements in strong electromagnetic interference environments or places with high vibration.
3. When weighing after not using the scale for a long time or after moving it, please step on the scale surface first, and wait for it to return to zero. Let it switch off automatically before switching on and weighing.
4. To ensure the accuracy of the fat measurement, please take off your shoes and socks, try to remove the influence of other weights outside your body, wash and wipe your hands and feet, and then step on the scale with your bare feet and hold the handle with both hands to take the measurement.
5. The scale surface is smooth. Please do not use this product on slippery ground or step on the edge of the scale. When weighing on the scale, first stand in the middle of the scale surface. Adjust your feet after standing steadily and step evenly on the 4 electrodes.
Please keep your body stable and balanced to prevent the side from tipping or slipping.
6. Do not measure after strenuous exercise, overeating, bathing, sauna, hot springs, or in a state of extreme dehydration.
7. The product should be kept away from high temperatures. Please do not use the product in humid, too hot, or too cold environment (using temperature range 5℃~35℃).
8. The product is a precision instrument, and cannot be disassembled. Please do not drop or knock.
9. Please keep the product clean. Do not use acidic detergent to clean it.

FAQ

The software can't obtain weighing results from scale.

- Please confirm the connection sign on the smart device to know the connection status.
If it's a WiFi scale, you may also need to check if the WiFi devices are switched on.
- Confirm the Bluetooth is on in your device and the scale is on within the valid range of signal coverage.
- Check whether the battery power is running low, please replace with a new battery.
- Relaunch the application.
- Turn off the Bluetooth and on again, then bridge the Bluetooth connection between scale and smart phone again.

Fail to connect to your device

- Please make sure only one device is setting up connection with the scale.
- Please make sure your device is equipped with Bluetooth 4.0 or above.

If the product fails to operate normally or displays any abnormalities:

- Step onto the scale to turn it on.
- If the display remains blank or the abnormalities persist, check the batteries and their charge level. If necessary, remove the batteries for 3 seconds to restart the product or replace them with new ones.
- If the issue persists after replacement, contact your local after-sales service team.

GUARANTEE

This article is guaranteed for 3 years, in the terms and conditions set out in Act 1/2007 (Spanish legislation).

For any queries related to the device purchased, consult El Corte Inglés, Hipercor or Supercor.

Assistance Tel. : (+34) 649 478 710 or (+34) 900 373 900.

Damage or faults caused by the following are excluded from the guarantee:

- Incorrect installation (voltage, gas or water pressure, electrical or water connections), reinstallations or housings made by the consumer without applying the correct instructions.
- Accidental causes such as falls, knocks, contact with liquids, insertion of foreign bodies, or any other cause of force majeure.
- Negligent, inadequate or non-domestic use, such as appliances installed in hairdressers, bars, restaurants, hotels, etc.
- Manipulation by technical services other than the official ones of the brand.
- Corrosion and/or rust caused by the normal wear and tear of the appliance or accelerated by adverse environmental conditions.
- Use of accessories or consumables which are not the original ones of the brand.

The following are also excluded from the guarantee:

- Components exposed to wear and tear due to normal use (light bulbs, seals, insulators, tubes, drains, etc.), from the sixth month, save manufacturing defects.
- Non-electromechanical components, aesthetic, plastic, glass or folding components, soap dishes, shelves, grilles, etc.
- Conservation services, cleaning, unblocking, change of direction of door, removal of foreign bodies, obstructions, recalibration or fine-tuning, etc.
- Computer products: Elimination of viruses, restoration of programs for this reason, or the reinstallation of the hard disk because it has been wiped.

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Hermosilla 112 - 28009 Madrid

España / Espanha / Spain. NIF: A-28017895

clientes@elcorteingles.es

Fabricado en / na China / Made in China

elcorteingles.es

clientes@elcorteingles.es

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El Corte Inglés S.A.