perifit

PERIFIT PEARLS

User Guide

support@perifit.co

1 English
4 Français
7 Deutsch
10 Italiano
13 Español
16 Português
19 Svenska
22 Nederlandse

Perifit - X6 Innovations
128 rue la Boetie, 75008 Paris, France
TF-IFU-Pearls-v20260529

EN

perifit

Perifit Pearls User Guide
English

1. Introduction
Welcome to the Perifit Pearls User Guide. Perifit Pearls, also known as Kegel Balls, are used to strengthen the pelvic floor muscles. This guide will help you use them safely and effectively.

2. Quick Start
2.1 Ensure Cleanliness:
Make sure the Perifit Pearls are clean and dry.
2.2 Apply Lubricant:
Use a water-based lubricant and gently insert the Perifit Pearls.
2.3 Perform Exercises:
Contract and relax your pelvic floor muscles to perform your Kegel exercises.
2.4 Remove Gently:
Relax and gently pull the retrieval string or bear down to remove.
2.5 Clean After Use:
Wash with warm water and mild soap; dry thoroughly before storing.

EN

3. Detailed Instructions

3.1 Preparation and Safety Tips

- Consult your healthcare provider before starting, especially if you have pelvic floor issues.
- Always clean the Perifit Pearls before and after use with warm water and mild soap.
- Use a water-based lubricant to facilitate easy insertion.

3.2 Insertion

- Find a comfortable position (lying down or standing with one leg raised).
- Apply a small amount of lubricant to the balls and your vaginal opening.
- Gently insert the first ball, followed by the second if using a dual set.

3.3 Exercises

- Start with basic Kegel exercises: contract your pelvic floor muscles as if stopping the flow of urine, hold for 5 seconds, then relax for 5 seconds. Repeat 10 times.
- Gradually increase the duration and intensity as your strength improves.

3.4 Removal

- Relax your body and gently pull the retrieval string or bear down to ease the balls out.

3.5 Cleaning and Maintenance

- Wash the balls with warm water and mild soap after each use.
- Dry thoroughly before storing in a clean, dry place.

4. Safety Warnings


- Do not use if you experience pain or discomfort.
- Not suitable for use during pregnancy without medical advice.
- Keep out of reach of children.

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